

Prevently.ai - Complete Technical Specification

Table of Contents

1. [Project Overview](#)
 2. [Technical Architecture](#)
 3. [Feature Specifications](#)
 4. [UI/UX Requirements](#)
 5. [Component Architecture](#)
 6. [Setup Instructions](#)
 7. [Critical Fixes Implemented](#)
 8. [API Integration](#)
 9. [Testing Requirements](#)
 10. [Deployment Specifications](#)
-

1. Project Overview

App Description and Vision

Prevently.ai is a comprehensive preventive healthcare platform that leverages AI to provide personalized health insights, medication management, and family health coordination. The app serves as a digital health twin, tracking user vitals, providing AI-generated health tips, and facilitating proactive health management.

Target Audience

- Health-conscious individuals aged 25-65
- Families seeking coordinated health management
- Patients with chronic conditions requiring medication tracking
- Users interested in preventive healthcare and wellness optimization

Core Value Proposition

- **AI-Powered Personalization:** Custom health tips based on user profile and health data
 - **Digital Health Twin:** Comprehensive health modeling and risk assessment
 - **Family Health Hub:** Coordinated health management for family members
 - **Preventive Focus:** Emphasis on prevention rather than treatment
 - **Professional UI/UX:** Liquid glass aesthetic with neumorphic design elements
-

2. Technical Architecture

Tech Stack with Versions

Frontend

- **Next.js:** 14.2.28 (App Router)
- **React:** 18.2.0
- **TypeScript:** 5.2.2
- **Tailwind CSS:** 3.3.3
- **Framer Motion:** 12.23.0 (animations)
- **Radix UI:** Various components (2.x.x)

Backend & Database

- **Prisma ORM:** 6.7.0
- **PostgreSQL:** Database
- **NextAuth.js:** 4.24.11 (authentication)
- **bcryptjs:** 2.4.3 (password hashing)

UI Components

- **Radix UI Primitives:** Complete component library
- **Lucide React:** 0.446.0 (icons)
- **React Hook Form:** 7.53.0 (form management)
- **Zod:** 3.23.8 (validation)

Charts & Visualization

- **Chart.js:** 4.4.9
- **React Chart.js 2:** 5.3.0
- **Plotly.js:** 2.35.3
- **React Plotly.js:** 2.6.0
- **Recharts:** 3.0.2

Additional Libraries

- **Date-fns:** 3.6.0 (date manipulation)
- **React Hot Toast:** 2.4.1 (notifications)
- **Zustand:** 5.0.3 (state management)
- **SWR:** 2.2.4 (data fetching)

Project Structure

app/		# Next.js App Router
app/		# API routes
api/		# Authentication endpoints
auth/		# Digital twin data
digital-twin/		# Family management
family/		# AI insights
insights/		# Medication tracking
medications/		# User onboarding
onboarding/		# Health surveys
surveys/		# Health tips
tips/		# User management
user/		# Dashboard page
dashboard/		# Digital twin visualization
digital-twin/		# Family management
family/		# Health knowledge base
knowledge-bank/		# Medication management
medications/		# User onboarding flow
onboarding/		# User settings
settings/		# Health tips browser
tips/		# Global styles
globals.css		# Root layout
layout.tsx		# Landing/auth page
page.tsx		# React components
components/		# Authentication components
auth/		# Dashboard-specific components
dashboard/		# Digital twin components
digital-twin/		# Family management components
family/		# Knowledge base components
knowledge-bank/		# Layout components
layout/		# Medication components
medications/		# Onboarding wizard
onboarding/		# Context providers
providers/		# Settings components
settings/		# Tips components
tips/		# Reusable UI components
ui/		# Utility libraries
lib/		# NextAuth configuration
auth.ts		# Prisma client
db.ts		# TypeScript types
types.ts		# Utility functions
utils.ts		# Database schema
prisma/		# Prisma schema
schema.prisma		# Database scripts
scripts/		# Database seeding
seed.ts		# Static assets
public/		# App logo
logo.svg		

Database Schema

Core Tables

Users Table

```

model User {
  id          String    @id @default(cuid())
  name        String?
  email       String    @unique
  emailVerified DateTime?
  image       String?
  password    String?

  // Onboarding fields
  onboardingCompleted Boolean @default(false)
  dateOfBirth         DateTime?
  gender              String?
  height              Float?
  weight              Float?
  activityLevel       String?
  healthGoals         String[]
  medicalConditions   String[]
  allergies            String[]
  currentMedications  String[]

  // Preferences
  notificationsEnabled Boolean @default(true)
  weeklyReportsEnabled Boolean @default(true)
  familySharingEnabled Boolean @default(false)
  aiPersonalizationLevel String @default("balanced")
  dataRetentionPreference String @default("standard")

  // Relations
  accounts Account[]
  sessions Session[]
  healthTips HealthTip[]
  medications Medication[]
  familyMembers FamilyMember[]
  digitalTwinData DigitalTwinData[]
  surveyResponses SurveyResponse[]
  accountabilityPartnerships AccountabilityPartnership[]

  createdAt DateTime @default(now())
  updatedAt DateTime @updatedAt
}

```

HealthTip Table

```

model HealthTip {
  id          String    @id @default(cuid())
  userId      String
  user        User       @relation(fields: [userId], references: [id])

  title       String
  content     String    @db.Text
  category    String    // nutrition, exercise, mental_health, sleep, hydration, general
  priority    String    @default("medium") // low, medium, high, urgent

  isPersonalized Boolean @default(false)
  viewed        Boolean @default(false)
  rating        Int?    // 1-5 star rating

  // AI and Expert Review
  aiGenerated   Boolean @default(false)
  expertReviewed Boolean @default(false)
  evidenceLevel String? // Research quality indicator
  sourceUrls    String[] // Links to research/sources

  createdAt    DateTime @default(now())
  updatedAt    DateTime @updatedAt
}

```

Medication Table

```

model Medication {
  id          String    @id @default(cuid())
  userId      String
  user        User       @relation(fields: [userId], references: [id])

  name        String
  dosage      String
  frequency   String    // daily, twice_daily, weekly, as_needed
  timeOfDay   String[]  // morning, afternoon, evening, night
  instructions String?
  prescribedBy String?

  startDate   DateTime
  endDate     DateTime?

  status      String    @default("active") // active, paused, completed, discontinued

  // Tracking
  adherenceRate Float?  // Percentage of doses taken as prescribed
  lastTaken    DateTime?
  nextDue      DateTime?

  createdAt    DateTime @default(now())
  updatedAt    DateTime @updatedAt
}

```

DigitalTwinData Table

```

model DigitalTwinData {
  id          String   @id @default(cuid())
  userId      String
  user        User     @relation(fields: [userId], references: [id])

  bodySystem  String   // cardiovascular, nervous, respiratory, digestive, immune, en-
doctrine
  dataType    String   // Heart Rate, Blood Pressure, Stress Level, etc.
  value       Float
  unit        String

  riskLevel   String   @default("low") // low, moderate, high
  trend       String   @default("stable") // improving, stable, declining

  // Data source
  source      String?  // manual, device, lab_result, app_sync
  deviceId    String?

  // Metadata
  notes       String?
  isVerified  Boolean  @default(false)

  lastUpdated DateTime @default(now())
  createdAt   DateTime @default(now())
}

```

API Endpoints

Authentication

- POST /api/auth/signup - User registration
- POST /api/auth/[...nextauth] - NextAuth.js handlers

Health Tips

- POST /api/tips/generate - Generate AI-powered health tip
- POST /api/tips/rate - Rate a health tip
- POST /api/tips/generate-category - Generate category-specific tip

Medications

- GET /api/medications - Get user medications
- POST /api/medications - Add new medication
- PUT /api/medications/[id] - Update medication
- DELETE /api/medications/[id] - Delete medication
- POST /api/medications/log - Log medication taken

Family Management

- GET /api/family - Get family members
- POST /api/family - Add family member
- PUT /api/family/[id] - Update family member
- DELETE /api/family/[id] - Delete family member
- POST /api/family/accountability - Create accountability partnership

Digital Twin

- GET /api/digital-twin/data - Get digital twin data
- POST /api/digital-twin/data - Add health data point

User Management

- GET /api/user/profile - Get user profile
- PUT /api/user/profile - Update user profile
- POST /api/user/water-intake - Log water intake
- GET /api/user/settings - Get user settings
- PUT /api/user/settings - Update user settings
- POST /api/user/export - Export user data
- DELETE /api/user/delete - Delete user account

Surveys

- POST /api/surveys/complete - Complete survey

Insights

- POST /api/insights/feedback - Provide feedback on AI insight

Onboarding

- POST /api/onboarding - Save onboarding data

3. Feature Specifications

Core Features

1. User Authentication

- **Demo Account:** john@doe.com / johndoe123
- **Secure Login:** bcrypt password hashing
- **Session Management:** JWT-based sessions via NextAuth.js
- **Registration Flow:** Email/password with validation

2. AI-Powered Health Tips

- **Personalized Generation:** Tips based on user profile, health data, and preferences
- **Categories:** Prevention, nutrition, exercise, mental health, sleep, hydration, general
- **Rating System:** 1-5 star rating for tip quality
- **Evidence-Based:** Tips include evidence level and source URLs
- **Fresh Content:** Daily generation of new personalized tips

3. Digital Health Twin

- **Body Systems Tracking:** Cardiovascular, respiratory, nervous, digestive, immune, endocrine
- **Risk Assessment:** Low, moderate, high risk levels
- **Trend Analysis:** Improving, stable, declining trends
- **Data Sources:** Manual entry, device integration, lab results
- **Visualization:** Interactive charts and health model

4. Medication Management

- **Comprehensive Tracking:** Name, dosage, frequency, timing
- **Adherence Monitoring:** Track missed doses and adherence rates
- **Smart Reminders:** Time-based medication alerts
- **Prescription Management:** Prescriber information and instructions
- **Status Tracking:** Active, paused, completed, discontinued

5. Family Health Hub

- **Family Member Profiles:** Name, relationship, age, health summary
- **Emergency Contacts:** Designated emergency contact system
- **Consent Management:** Privacy and data sharing controls
- **Accountability Partners:** Mutual health goal support
- **Sharing Levels:** Basic, detailed, full health data sharing

6. Interactive Health Surveys

- **Flo-Style Engagement:** Modern, engaging survey interface
- **Multiple Question Types:** Scale, multiple choice, text, boolean, slider
- **Daily Check-ins:** Mood, energy, sleep quality assessments
- **AI Insights:** Generated insights from survey responses
- **Progress Tracking:** Historical survey data and trends

7. Knowledge Bank

- **Clinical Studies:** Evidence-based health information
- **Expert Content:** Professionally reviewed health articles
- **Search Functionality:** Find relevant health information
- **Categorized Content:** Organized by health topics

8. Settings & Preferences

- **Dark Mode:** System-wide dark theme support
- **High Contrast:** Accessibility mode for better visibility
- **Notification Settings:** Granular notification controls
- **Privacy Controls:** Data retention and sharing preferences
- **AI Personalization:** Adjustable AI recommendation levels

User Flows

Onboarding Flow

1. **Personal Information:** Age, gender, height, weight, emergency contact
2. **Medical History:** Conditions, allergies, surgeries, family history
3. **Lifestyle Assessment:** Activity level, sleep, diet, stress
4. **Vital Signs:** Heart rate, blood pressure, BMI
5. **Health Goals:** Primary goals, timeline, motivation

Daily Usage Flow

1. **Dashboard Landing:** Health score, streak, today's tip
2. **Tip Interaction:** Read, rate, generate new tips
3. **Health Data Entry:** Log vitals, symptoms, mood
4. **Medication Tracking:** Mark medications as taken
5. **Survey Completion:** Daily health check-ins

Family Management Flow

1. **Add Family Member:** Basic information and relationship
 2. **Set Permissions:** Consent and sharing level configuration
 3. **Create Accountability:** Set up mutual health reminders
 4. **Monitor Progress:** View family health summaries
-

4. UI/UX Requirements

Design System

Liquid Glass Aesthetic

The app implements Apple's liquid glass design language with the following characteristics:

Core Glass Effects

```
.liquid-glass {
  background: rgba(255, 255, 255, 0.85);
  backdrop-filter: blur(20px) saturate(1.8);
  -webkit-backdrop-filter: blur(20px) saturate(1.8);
  border: 1px solid rgba(255, 255, 255, 0.3);
  box-shadow:
    0 8px 32px rgba(0, 0, 0, 0.1),
    inset 0 1px 0 rgba(255, 255, 255, 0.8),
    inset 0 -1px 0 rgba(255, 255, 255, 0.2);
}
```

Interactive Elements

```
.liquid-button {
  background: rgba(255, 255, 255, 0.9);
  backdrop-filter: blur(12px) saturate(1.4);
  border-radius: 12px;
  transition: all 0.2s ease;
}

.liquid-button:hover {
  transform: translateY(-2px);
  box-shadow: 0 8px 40px rgba(0, 0, 0, 0.12);
}
```

Neumorphic Design Elements

- **Soft Shadows:** Subtle depth with multiple shadow layers
- **Rounded Corners:** 12-16px border radius for cards and buttons
- **Pressed States:** Inset shadows for active elements
- **Gradient Backgrounds:** Subtle gradients for depth

Color Scheme

Light Mode

```
:root {
  --background: 0 0% 100%;
  --foreground: 240 10% 3.9%;
  --primary: 221.2 83.2% 53.3%;
  --secondary: 210 40% 96%;
  --muted: 210 40% 96%;
  --accent: 210 40% 96%;
  --border: 214.3 31.8% 91.4%;
}
```

Dark Mode

```
.dark {
  --background: 240 10% 3.9%;
  --foreground: 0 0% 98%;
  --primary: 221.2 83.2% 53.3%;
  --secondary: 240 3.7% 15.9%;
  --muted: 240 3.7% 15.9%;
  --accent: 240 3.7% 15.9%;
  --border: 240 3.7% 15.9%;
}
```

High Contrast Mode

```
.high-contrast {
  --background: 0 0% 0%;
  --foreground: 0 0% 100%;
  --primary: 0 0% 100%;
  --muted: 0 0% 20%;
}
```

Typography

- **Font Family:** System fonts (San Francisco, Segoe UI, Roboto)
- **Font Weights:** 400 (regular), 500 (medium), 600 (semibold), 700 (bold)
- **Font Sizes:**
 - Headings: 2xl (24px), xl (20px), lg (18px)
 - Body: base (16px), sm (14px)
 - Captions: xs (12px)

Responsive Design

- **Mobile First:** Designed for mobile with desktop enhancements
- **Breakpoints:**
 - sm: 640px
 - md: 768px
 - lg: 1024px
 - xl: 1280px
- **Grid System:** CSS Grid and Flexbox for layouts
- **Touch Targets:** Minimum 44px for interactive elements

Animation System


- **Framer Motion:** Page transitions and component animations
- **Micro-interactions:** Hover states, button presses, loading states
- **Breathing Animation:** Subtle scale animations for emphasis
- **Shimmer Effects:** Loading state animations
- **Spring Physics:** Natural feeling animations

Accessibility Features

- **High Contrast Mode:** Enhanced visibility for users with visual impairments
- **Focus Management:** Clear focus indicators and logical tab order
- **Screen Reader Support:** Proper ARIA labels and semantic HTML
- **Keyboard Navigation:** Full keyboard accessibility
- **Color Contrast:** WCAG AA compliant color ratios

5. Component Architecture

Component Tree



Reusable Component Specifications

NeumorphicCard

```

interface NeumorphicCardProps {
  children: React.ReactNode;
  className?: string;
  onClick?: () => void;
  variant?: 'default' | 'pressed' | 'elevated';
}

```

NeumorphicButton

```
interface NeumorphicButtonProps {
  children: React.ReactNode;
  onClick?: () => void;
  variant?: 'default' | 'primary' | 'secondary' | 'ghost';
  size?: 'sm' | 'md' | 'lg';
  disabled?: boolean;
  loading?: boolean;
}
```

SurveyComponent

```
interface SurveyComponentProps {
  survey: Survey;
  onComplete: (answers: Record<string, any>) => void;
  onClose: () => void;
}
```

AllInsightCard

```
interface AIInsightCardProps {
  insight: AIInsight;
  onFeedback: (insightId: string, helpful: boolean) => void;
  className?: string;
}
```

State Management Patterns

Local State

- **useState**: Component-level state for UI interactions
- **useEffect**: Side effects and data fetching
- **useReducer**: Complex state logic (survey responses, onboarding)

Global State

- **Zustand**: Lightweight state management for user preferences
- **SWR**: Server state management and caching
- **React Context**: Theme and authentication state

Form Management

- **React Hook Form**: Form state and validation
- **Zod**: Schema validation
- **Formik**: Complex multi-step forms (onboarding)

6. Setup Instructions

Prerequisites

- Node.js 18.x or higher
- PostgreSQL database
- Yarn or npm package manager

Environment Variables

Create `.env.local` file:

```
# Database
DATABASE_URL="postgresql://username:password@localhost:5432/prevently_db"

# NextAuth.js
NEXTAUTH_URL="http://localhost:3000"
NEXTAUTH_SECRET="your-secret-key"

# AI Service
ABACUSAI_API_KEY="your-abacus-ai-api-key"

# Optional: External integrations
GOOGLE_CLIENT_ID="your-google-client-id"
GOOGLE_CLIENT_SECRET="your-google-client-secret"
```

Step-by-Step Setup

1. Clone and Install Dependencies

```
# Clone the repository
git clone <repository-url>
cd prevently-ai/app

# Install dependencies
yarn install
# or
npm install
```

2. Database Setup

```
# Generate Prisma client
npx prisma generate

# Run database migrations
npx prisma db push

# Seed the database with demo data
npx prisma db seed
```

3. Development Server

```
# Start development server
yarn dev
# or
npm run dev
```

4. Build for Production

```
# Build the application
yarn build
# or
npm run build

# Start production server
yarn start
# or
npm start
```

Database Seeding

The seed script creates:

- Demo user: john@doe.com / johndoe123
- Sample medications (Vitamin D3, Omega-3)
- Health tips across all categories
- Family members with different relationships
- Digital twin data for all body systems
- Sample surveys and responses
- Accountability partnerships

Development Tools

- **ESLint:** Code linting with Next.js configuration
- **TypeScript:** Strict type checking enabled
- **Prettier:** Code formatting (via ESLint plugin)
- **Prisma Studio:** Database GUI (`npx prisma studio`)

7. Critical Fixes Implemented

Authentication Issues

Problem: Black screen on login, session management failures

Solution:

- Implemented proper NextAuth.js configuration with JWT strategy
- Added session callbacks for user ID persistence
- Fixed credential provider validation
- Added proper error handling for auth failures

UI Visibility Problems

Problem: Components not rendering, styling conflicts

Solution:

- Fixed CSS class conflicts between Tailwind and custom styles
- Implemented proper z-index layering for modals and overlays
- Added fallback styles for unsupported backdrop-filter
- Fixed responsive design breakpoints

TypeScript Errors

Problem: Type mismatches, missing interfaces

Solution:

- Created comprehensive type definitions in `lib/types.ts`
- Fixed Prisma client type imports
- Added proper prop interfaces for all components
- Implemented strict TypeScript configuration

Performance Optimizations

Problem: Slow page loads, unnecessary re-renders

Solution:

- Implemented SWR for efficient data fetching and caching
- Added React.memo for expensive components
- Optimized image loading with Next.js Image component
- Implemented lazy loading for non-critical components

Database Connection Issues

Problem: Prisma client connection failures

Solution:

- Fixed Prisma schema configuration for production
- Added proper connection pooling
- Implemented database connection retry logic
- Added environment-specific database URLs

Mobile Responsiveness

Problem: Poor mobile experience, touch targets too small

Solution:

- Implemented mobile-first responsive design
- Increased touch target sizes to 44px minimum
- Added proper viewport meta tags
- Optimized animations for mobile performance

Accessibility Improvements

Problem: Poor screen reader support, keyboard navigation issues

Solution:

- Added proper ARIA labels and roles
 - Implemented logical tab order
 - Added high contrast mode support
 - Fixed color contrast ratios for WCAG compliance
-

8. API Integration

AI Tips Generation System

Abacus.AI Integration

```
const response = await fetch('https://apps.abacus.ai/v1/chat/completions', {
  method: 'POST',
  headers: {
    'Content-Type': 'application/json',
    'Authorization': `Bearer ${process.env.ABACUSAI_API_KEY}`,
  },
  body: JSON.stringify({
    model: 'gpt-4.1-mini',
    messages: [{ role: 'user', content: prompt }],
    response_format: { type: "json_object" },
    temperature: 0.7,
    max_tokens: 500,
  }),
});
```

Personalization Algorithm

Tips are personalized based on:

- User age, gender, health conditions
- Current medications and allergies
- Lifestyle factors (activity level, sleep, stress)
- Previous tip ratings and preferences
- Digital twin health data

Content Quality Assurance

- Evidence-based recommendations only
- Source URL validation
- Expert review flags
- User rating feedback loop
- Content moderation for safety

External Service Integrations

Health Data Sources

- **Manual Entry:** User-input health data
- **Device Integration:** Fitness trackers, smart scales
- **Lab Results:** Import from healthcare providers
- **App Sync:** Integration with health apps

Notification Services

- **Email:** Weekly health reports
- **Push Notifications:** Medication reminders
- **SMS:** Emergency contact alerts
- **In-App:** Real-time health insights

Authentication Flow

NextAuth.js Configuration

```
export const authOptions: NextAuthOptions = {
  providers: [
    CredentialsProvider({
      name: 'credentials',
      credentials: {
        email: { label: 'Email', type: 'email' },
        password: { label: 'Password', type: 'password' },
      },
      async authorize(credentials) {
        // Validation and authentication logic
      },
    }),
  ],
  session: { strategy: 'jwt' },
  pages: { signIn: '/' },
  callbacks: {
    async jwt({ token, user }) {
      if (user) token.id = user.id;
      return token;
    },
    async session({ session, token }) {
      if (session.user) session.user.id = token.id as string;
      return session;
    },
  },
};
```

Data Validation and Error Handling

Input Validation

- **Zod Schemas:** Type-safe validation for all inputs
- **Sanitization:** XSS prevention and data cleaning
- **Rate Limiting:** API endpoint protection
- **CSRF Protection:** Built-in Next.js CSRF protection

Error Handling

```
try {
  // API operation
} catch (error) {
  console.error('Operation failed:', error);
  return NextResponse.json(
    { error: 'Operation failed' },
    { status: 500 }
  );
}
```

9. Testing Requirements

Unit Testing

- **Jest:** Testing framework for utility functions

- **React Testing Library:** Component testing
- **Coverage Target:** 80% code coverage minimum

Integration Testing

- **API Route Testing:** Test all endpoint functionality
- **Database Testing:** Prisma operations validation
- **Authentication Testing:** Login/logout flows

End-to-End Testing

- **Playwright:** Full user journey testing
- **Critical Paths:**
 - User registration and onboarding
 - Daily tip generation and rating
 - Medication tracking workflow
 - Family member management
 - Survey completion flow

Performance Testing

- **Lighthouse:** Performance, accessibility, SEO scores
- **Load Testing:** API endpoint stress testing
- **Mobile Performance:** Touch responsiveness testing

User Acceptance Criteria

Dashboard Functionality

- [] User can view personalized health dashboard
- [] Daily tip displays with rating capability
- [] Health score and streak are visible
- [] Water intake tracking works correctly
- [] Medication reminders appear when due

Onboarding Process

- [] New user can complete 5-step onboarding
- [] All form validation works correctly
- [] User profile is created successfully
- [] Initial health data is saved

Family Management

- [] User can add family members
- [] Consent and sharing levels work
- [] Accountability partnerships function
- [] Emergency contact system works

AI Tips System

- [] Tips generate based on user profile
 - [] Rating system functions correctly
 - [] Fresh tips appear daily
 - [] Evidence sources are provided
-

10. Deployment Specifications

Production Environment Setup

Hosting Requirements

- **Platform:** Vercel, Netlify, or AWS
- **Node.js:** 18.x or higher
- **Database:** PostgreSQL (managed service recommended)
- **CDN:** For static asset delivery

Environment Configuration

```
# Production Environment Variables
NODE_ENV=production
DATABASE_URL="postgresql://prod_user:password@prod_host:5432/prevently_prod"
NEXTAUTH_URL="https://prevently.ai"
NEXTAUTH_SECRET="production-secret-key"
ABACUSAI_API_KEY="production-api-key"
```

Build Configuration

```
// next.config.js
const nextConfig = {
  output: 'standalone',
  images: { unoptimized: true },
  eslint: { ignoreDuringBuilds: false },
  typescript: { ignoreBuildErrors: false },
};
```

Database Migration Strategy

1. **Backup:** Create database backup before deployment
2. **Migration:** Run `prisma db push` or `prisma migrate deploy`
3. **Seeding:** Optional production data seeding
4. **Verification:** Test database connectivity and queries

Monitoring and Logging

Application Monitoring

- **Error Tracking:** Sentry or similar service
- **Performance Monitoring:** Real user monitoring
- **Uptime Monitoring:** Health check endpoints
- **Database Monitoring:** Query performance tracking

Logging Strategy

- **Structured Logging:** JSON format for log aggregation
- **Log Levels:** Error, warn, info, debug
- **Sensitive Data:** Exclude PII from logs
- **Retention:** 30-day log retention policy

Security Considerations

Data Protection

- **Encryption:** All sensitive data encrypted at rest
- **HTTPS:** SSL/TLS for all communications
- **API Security:** Rate limiting and authentication
- **GDPR Compliance:** Data export and deletion capabilities

Privacy Controls

- **Data Minimization:** Collect only necessary data
- **Consent Management:** Granular privacy controls
- **Data Retention:** Configurable retention periods
- **Anonymization:** Option to anonymize health data

Backup and Recovery

Database Backups

- **Frequency:** Daily automated backups
- **Retention:** 30-day backup retention
- **Testing:** Monthly backup restoration tests
- **Geographic Distribution:** Multi-region backup storage

Disaster Recovery

- **RTO:** 4-hour recovery time objective
- **RPO:** 1-hour recovery point objective
- **Failover:** Automated failover procedures
- **Communication:** Incident response plan

Conclusion

This comprehensive specification provides everything needed to rebuild the Prevently.ai app in Cursor IDE. The application represents a modern, AI-powered healthcare platform with a focus on prevention, family coordination, and user engagement through beautiful UI/UX design.

Key implementation highlights:

- **Modern Tech Stack:** Next.js 14, React 18, TypeScript, Prisma
- **AI Integration:** Personalized health tips via Abacus.AI
- **Professional Design:** Liquid glass aesthetic with neumorphic elements
- **Comprehensive Features:** Digital twin, medication tracking, family hub
- **Production Ready:** Authentication, security, accessibility, performance

The app is currently functional and launch-ready, with all critical issues resolved and a complete feature set implemented.